

IMPORTANT WARNING WHEN REPLACING A CONTROLLER ON OPxxx, OMAXI, OMINI

On some old osmosers, the cable bundle has no white wire on the 4th pin of the connector.



old cable bundle:

standard OPxxx, OMAXI, OMINI cable bundle:



ANY CABLE BUNDLE WITHOUT WHITE WIRE ON PIN 4 MUST BE REPLACED BEFORE CONNECTING A NEW CONTROLLER !!!

(in such a case, the internal fuse will burn and the osmoser will not work)

CONNECTION OF THE ELECTRIC CABLE BUNDLE

[Original file](#)



CONNECTION OF THE ELECTRIC CABLE BUNDLE

WIRE	CONNECTED TO
------	--------------

Cable 1 Yellow	<input type="checkbox"/> Low Pressure Switch
Cable 2 Blue	<input type="checkbox"/> Cable connector of the High Pressure Switch
Cable 3 Red	<input type="checkbox"/> Low Solenoid Valve
Cable 4 Green Black	<input type="checkbox"/> <input type="checkbox"/> T jonction <input type="checkbox"/> <input type="checkbox"/>
Cable 5 Green White	<input type="checkbox"/> High Solenoid Valve
Cable 6 Blue Brown	<input type="checkbox"/> 24V Power supply

From: <https://diasys-technologies.com/waterpurifiers/> - **DiaSys Water Purifier Systems**

Permanent link: https://diasys-technologies.com/waterpurifiers/doku.php?id=replace_controllerwithoutwhitewire&rev=1573828084

Last update: **2019/11/15 14:28**