

**IMPORTANT WARNING WHEN REPLACING A CONTROLLER ON OPxxx, OMAXI, OMINI**

On some old osmosers, the cable bundle has no white wire on the 4th pin of the connector.



old cable bundle:

standard OPxxx, OMAXI, OMINI cable bundle:



---

**ANY CABLE BUNDLE WITHOUT WHITE WIRE ON PIN 4 MUST BE REPLACED BEFORE CONNECTING A NEW CONTROLLER !!!**



( in such a case, the internal fuse will burn and the osmoser will not work )

---

## **CONNECTION OF THE ELECTRIC CABLE BUNDLE**

[Original file](#)



WIRE	CONNECTED TO
<b>Cable 1 Yellow</b>	 <p data-bbox="432 1211 1037 1249">Low Pressure Switch</p>
<b>Cable 2 Blue</b>	 <p data-bbox="432 1850 1037 1888">Cable connector of the High Pressure Switch</p>

**Cable 3 Red**




Low Solenoid Valve

**Cable 4 Green Black**



T jonction

**Cable 5 Green White**



High Solenoid Valve

**Cable 6 Blue Brown**



24V Power supply

From: <https://diasys-technologies.com/waterpurifiers/> - **DiaSys Water Purifier Systems**

Permanent link: [https://diasys-technologies.com/waterpurifiers/doku.php?id=replace\\_controllerwithoutwhitewire&rev=1573833428](https://diasys-technologies.com/waterpurifiers/doku.php?id=replace_controllerwithoutwhitewire&rev=1573833428)

Last update: **2019/11/15 15:57**