2025/11/04 16:22

## **OSMOSERS-TROUBLESHOOTING**

## PROBLEMS POSSIBLES CAUSES SUGGESTED ACTION

- indicator *PUMP* switched off| - Switch on OFF position.| Switch ON the device|

	<ul> <li>The water supply of the tap closed or insufficient pressure (in that case the light "SOURCE" is on).</li> </ul>	Check the water supply
	- water hose folded or obstructed	Check that the inlet pipe is not damaged.
<ul> <li>No production of purified water - Power light ON - all control light indicators OFF</li> </ul>	- on power supply	Check and replace the fuses
Not an law assisted and a	Pre-treatment clogged	Replace the pre-treatment filters before the RO membrane
- Not or low purified water production - indicator <i>PUMP</i> light ON - Pump running	Low flow of inlet water	Increase the flow of the inlet water
light ON Tump running	Reverse osmosis membrane clogged	Change the membrane
- Not or low purified water production - indicator <i>PUMP</i> light ON - Pump always OFF		Replace the pump fuse
Light indicators FULL and PUMP alternate	Defectiveness of the high-pressure sensor or check valve.	Pressure sensor and check valve
	- Watertank hand valve closed.	Turn ON the watertank hand valve
Lack of pressure and low flow output	<ul> <li>volume of pure water need above the specification of the water purifier</li> </ul>	Wait the watertank is filled again
	<ul> <li>Lack of pressure in the pressurised watertank</li> </ul>	Re-adjust the air-pressure of the water tank
purified water conductivity too high	Resin saturated.	Replace resin cartridge or bottle
- POWER & FULL indicators ON - Pump OFF - Drain flowing continuously	- Inlet solenoid valve defective	Replace the inlet solenoid valve
Conductivity meter display OFF	Batteries out of order	Replace the batteries

## From:

https://diasys-technologies.com/waterpurifiers/ - DiaSys Water Purifier Systems

Permanent link:

https://diasys-technologies.com/waterpurifiers/doku.php?id=rotroubleshooting&rev=1554382855

Last update: 2019/04/04 13:00