Membrane replacement test

- Measure the conductivity after the reverse osmosis membrane
- Measure the conductivity of the tap water
- Do the following calculation: 100 * (conductivity reverse osmosis) / (conductivity raw water)
- If the value is less than 10%, the membrane is still usable.

In the other case, the reverse osmosis membrane begins to saturate.

From:

https://diasys-technologies.com/waterpurifiers/ - DiaSys Water Purifier Systems

Permanent link:

https://diasys-technologies.com/waterpurifiers/doku.php?id=when replacing membrane & rev=15536815321 and replaced to the control of the con

Last update: 2019/03/27 10:12