

# Conductivity Meter with Automatic Temperature Compensation

Reference DiaSys: **950186**

This easy to use and robust conductivity meter measures with accuracy the tap water as well as the deionised water produced by the water purifiers. It allows to control the system.



## SPECIFICATIONS



- Pen-type
- Large screen with backlight
- 4-digits LCD
- Automatic Temperature Compensation
- Measuring range: 0~1999us/cm



- Measurement accuracy:  $\pm 2\%$
- Power supply voltage: 4x1.5v (AG-13 button battery)
- Resolution : 1us/cm



{{:wiki:ec2013-3.jpg?600 |}}  
{{:wiki:ec2013-4.jpg?600 |}}

From:  
<https://diasys-technologies.com/waterpurifiers/> - **DiaSys Water Purifier Systems**

Permanent link:  
<https://diasys-technologies.com/waterpurifiers/doku.php?id=ec2013&rev=1553528616>



Last update: **2019/03/25 15:43**