

Range of water purifiers

O-MINI

- up to 30 liters per day
 - water tank of 12 liters
 - Volume of resin filter : 0.75 liter
 - Especially convenient for a clinical analyser up to 200T/H
-

O-MAXI

- up to 50 liters per day
 - water tank of 40 liters
 - Volume of resin filter : 12 liters
 - Especially convenient for a clinical analyser up to 600T/H
-

user manuals

OMINI user manual version February 2019

OMAXI user manual version February 2019

From:

<https://diasys-technologies.com/waterpurifiers/> - **DiaSys Water Purifier Systems**

Permanent link:

<https://diasys-technologies.com/waterpurifiers/doku.php?id=waterpurifiers&rev=1553682616> 

Last update: **2019/03/27 10:30**