

# Membrane replacement test

Do the following calculation to get a percentage:  $100 * (\text{TDS reverse osmosis}) / (\text{TDS raw water}) = x\%$

If the value is less **than 10%**, the membrane is still usable. In the other case, the reverse osmosis membrane begins to saturate.

From:

<https://diasys-technologies.com/waterpurifiers/> - **DiaSys Water Purifier Systems**

Permanent link:

<https://diasys-technologies.com/waterpurifiers/doku.php?id=whenreplacingmembrane&rev=1553681404> 

Last update: **2019/03/27 10:10**